

Inhaler Technique

1. Open the case until it clicks
2. Slide lever back as far as it can go until it clicks
3. Breathe out fully
4. Seal your lips around the mouthpiece
5. Breathe in slowly, steadily and deeply
6. Hold your breath for about 10 seconds or for as long as is comfortable
7. Remove the device from your mouth
8. Close the case
9. Repeat steps 1-7 for each dose of medication prescribed

Cleaning instructions:

Never wash inhaler with water, just wipe the mouthpiece and the outside with a clean, dry, lint-free cloth to remove any powder residue.

Keep inhaler dry.

Further advice:

Always remember to gargle with water after using your inhalers and spit out the residue.

You can request this information in a different format or another language.

Date of publication: July 2013,
reviewed Oct 2016

Author: Clare Walsh CNS, Claire Wotton CNS

Reference Number: ACC v1


Review Date: Oct 2018

© East and North Hertfordshire NHS Trust

www.enherts-tr.nhs.uk

Department of Respiratory Medicine

Accuhaler information leaflet



For further information contact the
Respiratory Specialist Nurses
Lister Hospital: **01438 285621**

