

## Inhaler Technique

1. Remove the protective cap by lightly squeezing the arrows
2. Press the green button all the way down
3. Release the green button
4. Check that the colour control window is green. It means the inhaler is ready
5. Inhale **STRONGLY** and **DEEPLY** through the mouthpiece. Keep breathing in, even after you have heard the inhaler "click"
6. Check that the colour control window has turned to red from green, which indicates you have inhaled correctly
7. Replace the dust cap.

### INCORRECT USE;

- \* **DO NOT BLOW INTO THE INHALER**
- \* **DO NOT HOLD THE GREEN BUTTON DOWN WHILE INHALING**
- \* **CONTROL WINDOW REMAINS GREEN AFTER INHALATION**

If the colour control window is still green after inhalation, please repeat inhaling **STRONGLY** and **DEEPLY** through the mouthpiece.

### Further advice:

This inhaler should be cleaned once a week. This can be achieved by:

Wipe the mouthpiece inside and outside with a clean, dry, lint-free cloth to remove any powder residue.

Never wash the inhaler with water, just wipe the mouthpiece

Keep inhaler dry.

You can request this information in a different format or another language.

Date of publication: July 2013,  
reviewed Oct 2016

Author: Clare Walsh CNS, Claire Wotton CNS

Reference Number: GEN v1

Review Date: Oct 2018

© East and North Hertfordshire NHS Trust

[www.enher-ts-tr.nhs.uk](http://www.enher-ts-tr.nhs.uk)

## Department of Respiratory Medicine

### Genuair Inhaler information leaflet



For further information contact the  
Respiratory Specialist Nurses  
Lister Hospital: **01438 285621**

